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TP N°4: SILVER DETERMINATION OF CHLORIDES

Precipitation methods are based on reactions involving poorly soluble compounds. In these methods, the equivalence point is characterized by the complete precipitation of the ions to be determined.

Although many reactions are known to be accompanied by the formation of insoluble compounds, only some of them can be used in the chemical industry. Only some of them can be used in volumetric analysis. These reactions must satisfy a number of conditions:

- The precipitate must be practically insoluble.
- Precipitation must occur with sufficient rapidity (the formation of a supersaturated solution is absent).
- It must be possible to set the equivalence point during titration.

The most important precipitation methods are :



1. Argentimetry: based on the following reaction : AgX with X^- : Cl^- , Br^- , I^-

2. Mercurimetry: based on the following reaction: $\text{Hg}^{++} + 2\text{X}^-$.

Principle of silver chloride assay Argentimetric

assays Argentimetric are based on the property of silver salts, when mixed with an antagonistic substance, to produce an insoluble compound that precipitates out of solution.

Two methods are used:

a. Gay's method : In this method, the reaction is complete when the addition of excess Ag^+ no longer produces a precipitate. This means that all the Cl has precipitated.


b. Mohr's method : This method consists in introducing an end-of-reaction indicator, potassium chromate (K_2CrO_4).

As silver chromate is more soluble than silver chloride, it only precipitates once all the silver has been added.

The reaction ends with the addition of the drop that precipitates a red precipitate of silver chromate (Ag_2CrO_4). is precipitated.


Procedure

a. Gay's method

1. Take 10ml of (NaCl) of unknown molarity into an Erlenmeyer flask.
2. Fill the burette with a 0.1N solution of silver nitrate (AgNO_3).
3. Perform the first titration. Stop when the addition of a drop of AgNO_3 does not produce a precipitate.
4. Determine the volume of AgNO_3 used.
5. Perform another titration.

b. Mohr's method

1. Take 10ml of (NaCl) in an Erlenmeyer flask.
2. Add 2 drops of a cold saturated solution of K_2CrO_4 .
3. Fill burette with 0.1N (AgNO_3) solution.
4. Perform first titration. Stop as soon as a red precipitate appears.
5. Determine the volume of AgNO_3 used.
6. Perform another titration.