TP N°3

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.....

PRACTICAL EXERCISE 3

Determining the specific heat of a metal sample

Definitions

Calorimetry : Calorimetry is an operation used to measure the quantities of heat exchanged during transformations undergone by a system, as well as the determination of heat capacities.

Calorimeter: The calorimeter is an adiabatic vessel used in various calorimetry operations.

Calorimeter water value: In calorimetry the calorimeter water value is the fictitious mass of water μ that has the same heat capacity as the calorimeter.

$$m_{\text{calorimeter}} \cdot C_{\text{calorimeter}} = \mu \cdot C_{\text{eau}}$$

The water value of the calorimeter is a characteristic of the calorimeter, it must be determined experimentally and will be used for all manipulations.

$$C_{\text{eau}} = 4185 \text{ J/Kg/K}$$


Description of the adiabatic calorimeter : An adiabatic calorimeter (e.g. Berthelot type) consists of :

- 1 Vase made of insulating material (e.g. Dewar* vase)
- 2 Thermometer
- 3 Stirrer

* **Dewar vase :** The vase consists of two nested vessels with glass insulating walls, separated by a vacuum.

The near-vacuum prevents heat transfer by conduction and convection.

The inner surface of the outer vessel and the outer surface of the inner vessel have a reflective metallic coating to prevent heat being transferred by radiation. Silver is most often used for this purpose.

Principle

A mass of water m_1 is introduced into the calorimeter at room temperature T_1 , the same temperature as the calorimeter.

Add a sample of mass m_2 and temperature T_2 , and wait for equilibrium and the final temperature.

We apply the first principle of thermodynamics to the system:

(Calorimeter + Mass of water + Sample)

The calorimeter is an adiabatic system (no heat exchange with the outside) $Q_{\text{total}} = 0$.

$$Q_{\text{calorimeter}} + Q_{\text{water}} + Q_{\text{sample}} = 0$$

With

$$Q_{\text{calorimeter}} = \mu \cdot C_{\text{eau}} (T_F - T_1)$$

$$Q_{\text{water}} = m_1 C_{\text{eau}} (T_F - T_1)$$

$$Q_{\text{sample}} = m_2 \cdot C_{\text{sample}} (T_F - T_2)$$

Procedure

Introduce a mass of water $m_1=180$ g approximately at room temperature T_1 into the calorimeter. The calorimeter itself is at temperature T_1 . Measure temperature at regular intervals.

- Immerse a piece of metal of mass m_2 (to be determined using the balance) in boiling water, and hold it suspended using tongs or pliers for at least 5 minutes, to acquire a temperature of 100°C . Avoid contact of the metal sample with the hot plate.
- Introduce the metal sample, heated to 100°C , into the calorimeter.
- Shake the mixture and take the temperature at regular intervals. Record measurement results in the table.

Purpose of the second manipulation:

.....

Mass of water

Value of water in calorimeter

Mass of metal part

$m_{\text{water}} = m_1 =$

$\mu =$

$m_{\text{metal}} = m_2 =$


Calorimetric balance (equation)

Write the balance equation after quoting the calorific contribution of each component

.....

Table-1: Temperature variation as a function of time.

t (s)														
T (°C)														

Temperature of calorimeter + water $T_c = T_e =$

Temperature of metal part $T_m =$

Final temperature T_F (graph) =

Write down the analytical expression for C_{metal} (Heat capacity of the metal part):

.....

From the T_F measurement, determine the value of the metal's mass heat capacity C_{metal} (in $J/Kg/^\circ K$)

.....