TP CHIMIE1

Manipulation N° 5

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
Oxidation-reduction titration

General :

Oxidation-reduction reactions are electron exchange reactions involving electron acceptors and donors: oxidizing agents can fix electrons, while reducing agents can give them up.


The reaction that defines the relationship between a reducing agent and an oxidizing agent is :


In these reactions, the **reducing agent** oxidizes by gaining electrons, while the **oxidizing agent** reduces by losing electrons.

- The passage of iron "Fe" in solution corresponds to **oxidation**. $\overset{\sim}{\text{Fe}} \longrightarrow \text{Fe}^{2+} + 2 e$


- The deposition of copper "Cu" from the Cu²⁺ ion solution corresponds to a **reduction**.

Oxidant: is the reagent capable of causing oxidation, i.e. it captures one or more electrons from another reagent.

Reductant: is the reagent capable of causing a reduction, i.e. giving up one or more electrons to another reagent.

Redox couple: A redox couple is a couple formed by an oxidant and its conjugate reductant.

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Determination of iron (II) ions by permanganate ions in an acid medium

1. Aim:

To determine the normality of a solution (FeSO₄), using the oxidation reaction of a solution of Potassium Permanganate (KMnO₄) prepared in the laboratory in an acid medium. This assay is called manganimetry.

TP CHIMIE1


- **Manganimetry**: This method is based on permanganate ion oxidation reactions. Oxidation can be carried out in acid, alkaline or neutral media. However, acidic media are more often used in volumetric analysis.

- **Potassium permanganate** : Potassium permanganate (KMnO_4) is a particularly powerful oxidizing agent. It takes the form of violet crystals composed of potassium ions, K^+ , and permanganate ions, $[\text{MnO}_4^-]$ It is odorless and bitter-tasting.

In the laboratory, potassium permanganate is used for titration. At equivalence, the solution changes color from violet to pink. Potassium permanganate is also used in water treatment to oxidize iron and manganese in groundwater. In everyday life, it can be used in dilute solution to remove black marks left by fungi between bathroom tiles.

2. Handling principle

Dosing consists in determining the normality of a reducing solution given the normality of the oxidizing solution. We propose to study the oxidation of the Fe^{2+} ion by the permanganate ion MnO_4^- in an acid medium.

The oxidizing properties of the permanganate ion are the basis of manganometry. The oxidizing form, MnO_4^- , is violet, while the reducing form, Mn^{2+} , is colorless, enabling the equivalent point to be determined without the use of color indicators. H^+ ions are added in excess. Sulfuric acid H_2SO_4 is used. Excess sulfuric acid avoids reaction in a neutral medium.

3. Experimental part:

1. Take 10 ml of FeSO_4 of unknown normality N into an Erlenmeyer flask, and add approximately 5 ml of 2N H_2SO_4 .
2. Fill burette with KMnO_4 oxidizing solution of normality 0.1N.
3. Drip the KMnO_4 oxidizing solution with stirring. The first drop of KMnO_4 solution in excess, i.e. no longer reduced, colors the solution in the Erlenmeyer flask an evergreen pink.
4. Carry out 3 dosing tests.